CareFusion

# MicroCO Meter

Operating Manual – English

Operating Manual – CareFusion USA

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## Contents

I.	Introduction – English.....	4
	Smoke Inhalation/ CO Poisoning .....	5
	Cigarette Consumption .....	6
	Package Contents.....	7
	Operation.....	9
	PC Connection .....	11
	Calibration .....	12
	Battery Life .....	15
	Cleaning .....	16
	Servicing.....	17
	Consumables/ Supporting products .....	17
	Specifications .....	19
II.	Introduction – CareFusion USA .....	22
	Smoke Inhalation/ CO Poisoning .....	23
	Cigarette Consumption .....	24
	Package Contents.....	25
	Operation.....	27
	PC Connection .....	29
	Calibration .....	30
	Battery Life .....	33
	Cleaning .....	35
	Servicing.....	35
	Consumables/ Supporting products .....	35
	Specifications .....	37
III.	Introduction – Français .....	39
	Tabagisme / Empoisonnement au CO.....	40
	Tabagisme .....	41
	Contenu de l’emballage.....	42
	Consignes d’utilisation .....	44

Connexion PC .....	46
Calibration .....	46
Durée de vie de la pile .....	50
Nettoyage .....	51
Entretien.....	52
Produits de consommation/ autres produits .....	52
Spécifications .....	54
 IV. Einführung – Deutsch .....	 56
Rauchinhalation/ CO Vergiftung .....	57
Zigarettenkonsum .....	58
Packungsinhalte .....	60
Bedienung.....	62
PC Anschluss .....	65
Eichung.....	65
Lebensdauer der Batterie .....	70
Reinigung.....	71
Wartung.....	71
Verbrauchs-/ Hilfsmaterial.....	72
Anforderungen .....	73

# MicroCO Meter

## Operating manual

### Overview

The MicroCO meter is a hand held battery operated device used to measure the concentration of carbon monoxide, CO, on the breath and calculates the percentage of carboxyhaemoglobin, %COHb, in the blood.

It is used in smoking cessation clinics, GP's surgeries, accident and emergency departments, and by the fire-fighting services. In smoking cessation clinics it is used for instructional purposes and to check on the client's progress and compliance. General practitioners, A&E clinicians and the fire-fighting services can use the instrument to quickly assess the level of suspected CO poisoning. It is accurate, easy to use, and has many features designed to simplify its operation.

These include:

- Auto zero function
- Smoking level colour light indicators
- Poison level alarm
- Simple calibration
- Serial interface to PC

**Please note: It is recommended that this unit be calibrated upon receipt**

## I. Introduction - English

The MicroCO meter is based on an electrochemical fuel cell sensor, which works through the reaction of carbon monoxide (CO) with an electrolyte at one electrode, and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to CO concentration. Output from the sensor is monitored by a microprocessor, which detects peak expired concentrations of alveolar gas. This is then converted to % carboxyhaemoglobin (%COHb) using the mathematical relationships described by Jarvis et al, for concentrations below 90 ppm, and by Stewart et al for higher levels. Raised levels of carboxyhaemoglobin are most commonly caused by accidental smoke inhalation/CO poisoning or cigarette smoking. The MicroCO includes a countdown timer as an aid to timing the breath holding period prior to exhalation results are displayed on a clear LCD display. Warning lights are provided to give an instant indication of the smoking level.

The countdown timer, warning light levels, and the alarm level are user adjustable when the unit is connected to a PC running COBRA software.

**Note: the countdown timer, warning light levels, and the alarm level used throughout this manual are the factory settings and may have been changed.**

## References

**Jarvis MJ, Belcher M, Vesey C, Hutchison DCS**

**Low cost carbon monoxide monitors in smoking assessment**

Thorax 1986; 41: 886-887

Stewart RD, Stewart RS, Stamm W, Seleen RP

**Rapid estimation of carboxyhaemoglobin levels in fire fighters**

JAMA 1976; 235, 390-392

## Smoke Inhalation/CO Poisoning

CO poisoning usually occurs as a result of smoke inhalation from fires or exposure to CO from car exhausts or faulty heating systems. CO binds with haemoglobin to form carboxyhaemoglobin (COHb) which reduces the capacity of the blood to carry oxygen. Acute poisoning may cause symptoms ranging from headache and breathlessness (at COHb levels of 10% to 30%) through confusion to coma and death (COHb usually greater than 60%). A particularly insidious feature of CO poisoning is the development of neurological problems such as movement disorders (often resembling Parkinson's disease), memory loss and altered personality. Such problems may develop weeks after apparent recovery from acute poisoning.

Chronic exposure to relatively low levels of CO may result in a variety of symptoms including headache, fatigue, poor concentration, dizziness, palpitations, chest pain, visual disturbance, nausea, diarrhoea and abdominal pain.

Chronic CO poisoning as a cause of such symptoms is often missed.

The MicroCO meter allows immediate assessment of patients at risk of CO poisoning who can then be rapidly referred for expert assessment. Prompt treatment with oxygen (in a hyperbaric chamber, if severe) is often life saving. Screening for CO exposure may also reveal the cause of non-specific symptoms relating to low level CO exposure (usually as the result of faulty gas appliances).

## Reference

Meredith T, Vale A, **Carbon monoxide poisoning**  
British Medical Journal, 1988; 296, 77-78

## Cigarette Consumption

The MicroCO meter provides a simple screening test for cigarette consumption for use in anti-smoking clinics and all smoking cessation programs.

Measurement of carboxyhaemoglobin has been well validated as an indirect measure of cigarette consumption and is widely used in smoking cessation programs.

Typical values for carboxyhaemoglobin and expired CO in smokers and poisoning victims, together with the alarm light status, are given below:

CO(ppm)	%COHb	Cigarette consumption	Indicator
0 – 6	0 – 1	Non smoker	Green
7 –10	1.1 – 1.6	Light smoker	Amber
11 – 72	1.8 – 12	Heavy smoker	Red
>72	>12	Suspected poisoning	Red + alarm

Please note that some urban areas may have high environmental levels of CO. This can cause a rise in exhaled CO of a few ppm above that which is normally present on the breath. In these cases it is possible for a non-smoker to appear at the bottom of the 'light smoker' range (7 – 10 ppm).

## References

### **The relationship between alveolar and blood carbon monoxide concentrations during breath holding**

Jones RH, Ellicott MF, Cadigan JB, Gaensler EA

Journal of Laboratory and Clinical Medicine 1958; 51, 553 – 564

## **Carbon monoxide in breath in relation to smoking and carboxyhaemoglobin levels**

Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369

## **Definition of a reliable threshold value for detecting current smokers by CO measurement**

Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano - Respiratory and TB Unit-General  
Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-ITALY).

## **Correlation between exhaled CO measurements and carboxyhaemoglobin percentage in smokers**


Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla; Zanette  
Antonia; Nardini Stefano - Respiratory and TB Unit- General  
Hospital- Via Forlanini, 71- I-31029-Vittorio Veneto (TV- ITALY).

## **Package Contents**


The MicroCO meter is supplied with a carrying case containing the following items:

1. MicroCO meter (Cat No. 36-MC02-STK)
2. 22 mm mouthpiece adapter (Box of 10 Cat No.36-PSA2000)
3. 22 mm reducing connector (Cat No.36-MEC1007)
4. PP3 Battery (Cat No.36-BAT1002)
5. 4 Disposable mouthpieces (Box 250 Cat No.36-PSA1200)
6. Calibration tool (Cat No. 36-MEC1184)


1

2

3

4

5

6


## Operation


For accurate results the CO meter should be used at room temperature. If the instrument has been stored in cool or hot conditions then allow time to reach room temperature prior to use

Install the PP3 battery by sliding open the battery cover, clipping the battery in place and replacing the cover. Insert the mouthpiece adapter into the MicroCO meter and then insert a disposable mouthpiece into the adapter. Turn the unit on by selecting the upper or middle position on the central slide switch and the display will show the software version number.

This will appear momentarily whilst correction for ambient levels is executed. During this time the unit must not be exposed to raised concentrations of CO. The buzzer will then sound and the display will change to:

Upon hearing the buzzer instruct the subject to inspire fully and hold their breath for 20 seconds. The display will count down from 20 to 0 as an aid to timing the breath holding period. The green indicator will then illuminate and the display will change to:

At this time the subject should seal their lips around the mouthpiece and exhale slowly and fully. The 20 second breath holding period is recommended to allow time for equilibration of alveolar gas.

If, however, the subject is unable to hold their breath for this period, the unit may still be used before 20 seconds have elapsed.

The unit must not be used, however, for one second after the unit is first switched on i.e. before the countdown commences. Expired alveolar gas is then entrapped between sensor and mouthpiece valve. The display reading will rise to a plateau over the course of a several seconds. The final value will be held until the unit is turned off and represents parts per million CO or %COHb depending on the slide switch position.

The lights will come on according to the table on page 4 and the red light will flash and an alarm will sound if the measurement rises above 72ppm (12%COHb).

If this occurs then the possibility of CO poisoning should be investigated as this level of CO is unlikely to be produced by cigarette smoking.

### **Important note:**

**Before repeating a measurement the unit must be turned off, and the mouthpiece and adapter removed for at least 1 minute. This is to allow re-equilibration with ambient air and to dry the surface of the sensor. Visually inspect that all moisture has evaporated from the surface of the sensor before reuse.**

If the unit is switched on again too quickly after use there may be a response to residual expired carbon monoxide from the previous test.

In this case the display will show:

If this is displayed then turn the unit off, remove the mouthpiece adapter, and expose to ambient air for 2 minutes before repeating the test.

**Note:** If this warning appears again after following the above procedure then turn off and leave the sensor exposed to ambient air for a further 3 minutes. If the same message appears again then this indicates possible contamination of the sensor with a solvent.

In this case remove all sources of solvent from the vicinity of the sensor and expose to ambient air for 24 hours before switching on again.

## PC Connection

The MicroCO may be connected to the serial port of a PC running COBRA (CO BReath Analysis) software.

This software allows the measurements to be read by a PC and automatically entered onto a pre-defined report for subsequent printing and filing. It also allows the CO level indicators and the breath holding countdown timer to be configured to individual requirements.


**Note:** the MicroCO should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 – 'Safety of

Information Technology Equipment including Electrical Business Equipment’.

## Calibration


Calibration will remain stable to within 2% over one month and typically to within 10% over 6 months. CareFusion supplies calibration gas (20ppm CO in air) and recommends that the unit is recalibrated on a 6 monthly basis. See page 17 for calibration accessories.

To carry out the calibration locate the calibration switch on the right hand side of the instrument as shown below.

Push the slide switch to the CO - PPM position and wait for the unit to display zero.

Screw the control valve firmly onto the cylinder and connect the gas supply as shown below:

The plastic tubing supplied with the gas should be pushed firmly over the reducing connector.

Slowly turn the control knob anti-clockwise until the ball in the flow indicator is between the two marks. This will then supply a gas flow

of approximately 0.25 l/min. Apply this flow for 25 seconds and if the meter does not read 20ppm use the calibration tool to push the calibration switch.

The unit will then beep 3 times, store the new calibration value, and display the following:

The gas supply should then be turned off.


If the signal from the CO sensor is too low a new calibration value will not be stored and the display will show:

The most likely cause for this is an expired cell but may also be caused by depressing the calibration switch with no gas applied. Ensure that the concentration of calibration gas is correct (20ppm), that the connections to the gas cylinder are secure, and that the gas cylinder is not empty, and then repeat the calibration procedure.

If the above message is repeated, the CO meter must be returned to CareFusion, or an authorised agent, for sensor replacement. Sensor life ranges from 2 to 5 years and depends upon both the amount of exposure to CO and other gases, particularly solvents such as alcohol and cleaning fluids.

If the signal from the CO sensor is too high, a new calibration value will not be stored and the following will be displayed when calibration is attempted:

The most likely cause for this is using an incorrect concentration of calibration gas. Ensure that the concentration of calibration gas is correct (20ppm) that the connections to the gas cylinder are secure, and then repeat the calibration procedure.

In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.


### **Important Notes:**

- Only certified calibration gas from a reputable source should be used.
- Ensure that no CO is present on the sensor for 3 minutes before starting the calibration procedure.
- Ensure that the instrument and gas cylinder have stabilised at room temperature before calibrating.
- Do not immerse the unit in the calibration gas.

## **Battery Life**

Battery life is approximately 30 hours of continuous use.

When the battery has approximately 1 hour of useful life left the alarm will sound momentarily after the unit is first switched on and the following message will be displayed:


When the battery has completely expired the above will be displayed continuously and the battery must be replaced.

**Note:** Please remove the battery if the meter is likely to be unused for more than a month.

### Internal Battery Expiry

The MicroCO meter has an internal battery with a life of approximately 10 years. This battery supplies the sensor signal conditioning circuit continuously to ensure instant start up.

When the battery has expired, the following warning message will be displayed:

and the alarm will sound when the unit is first turned on.

When this occurs, the CO meter must be returned to CareFusion, or an authorised agent, for battery replacement.

### Cleaning

The mouthpiece adapter may be cleaned using a mild detergent solution or replaced if there is risk of contamination. Exposed

surfaces of the meter, other than the sensor surface, may be wiped with a clean, damp cloth. It is recommended that this procedure is performed after each use and all used cardboard mouthpieces discarded.

**Important note:** The sensor surface must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol or permanent damage may result.

## Servicing

If your unit requires service or repair please see page 75 for contact details. A full service manual including circuit diagram and parts list is available upon request.

## Consumables / Supporting Products

Cat. No.	Description
36-PSA1200	22mm disposable cardboard mouthpieces (250 per box)
36-PSA2000	22 mm mouthpiece adaptors (pack of 10)
36-BAT1002	Alkaline PP3 Battery
36-MEC1184	Calibration Tool
36-MCG020	MicroCan calibration gas. 20 litres of gas in a 1 litre canister. 20ppm Carbon Monoxide, balance air.
36-MGA222	Stainless steel valve for calibration canisters with flow Indicator
36-MEC1007	22mm reducing connector for calibration
36-CAB1000	Interface cable for COBRA software

## Caution

- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.”

## Symbols

Type B device

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

## Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

## Specifications

Sensor type	Electro-chemical fuel cell
Range	0 - 100 ppm
Resolution	1 ppm
Green indicator light	0 to 6ppm (0 to 1 %COHb)
Amber indicator light	7 to 10ppm (1.1 to 1.6 %COHb)
Red indicator light	11 to 72ppm (1.8 to 12 %COHb)
Flashing red light + alarm	>72ppm (>12 %COHb)
Accuracy	+/-5% of reading or 1ppm whichever is the greater
Sensitivity drift	0.5%/°C
Sensor life	2 to 5 years
Response time	< 15 sec (to 90% of reading)
Hydrogen cross sensitivity	<15%
Operating temperature	15 – 25 °C
Operating pressure	Atmospheric +/- 10%
Pressure coefficient	0.02% signal per mBar
Relative humidity (Non condensing)	15 - 90% continuous (0 - 99% intermittent)
Baseline drift	0ppm (auto-zero)
Long term drift	< 2% signal loss per month
Power source	Single Alkaline 9 volt PP3
Main battery life	30 hours of continuous use equivalent to approximately 2000 tests
Internal battery life	10 years
Weight	180 g (Including battery)
Dimensions	170 x 60 x 26 mm
Display	3 ½ digit LCD
Storage temperature	-20° to +70 °C
Storage humidity	30% to 90% RH

CareFusion

# MicroCO Meter

## Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,  
22745 Savi Ranch Parkway,  
Yorba Linda  
CA 92887-4668  
USA

## Overview

The MicroCO meter is a hand held battery operated device used to measure the concentration of carbon monoxide, CO, on the breath and calculates the percentage of carboxyhemoglobin, %COHb, in the blood.

It is used in smoking cessation clinics, physician offices, emergency departments, and by the fire-fighting services.

In smoking cessation clinics it is used for instructional purposes and to check on the client's progress and compliance. General practitioners, Emergency room personnel and fire-fighting services can use the instrument to quickly assess the level of suspected CO poisoning. It is accurate, easy to use, and has many features designed to simplify its operation.

These include:

- Auto zero function
- Smoking level color light indicators
- Poison level alarm
- Simple calibration
- Serial interface to PC

## II. Introduction – CareFusion USA

The MicroCO meter is based on an electrochemical fuel cell sensor, which works through the reaction of carbon monoxide (CO) with an electrolyte at one electrode, and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to CO concentration. Output from the sensor is monitored by a microprocessor, which detects peak expired concentrations of alveolar gas. This is then converted to % carboxyhemoglobin (%COHb) using the mathematical relationships described by Jarvis et al, for concentrations below 90 ppm, and by Stewart et al for higher levels. Raised levels of carboxyhemoglobin are most commonly caused by accidental smoke inhalation/CO poisoning or cigarette smoking.

The MicroCO includes a countdown timer as an aid to timing the breath holding period prior to exhalation.

The results are displayed on a clear LCD display. Warning lights are provided to give an instant indication of the smoking level.

The countdown timer, warning light levels, and the alarm level are user adjustable when the unit is connected to a PC running COBRA software.

**Note: the countdown timer, warning light levels, and the alarm level used throughout this manual are the factory settings and may have been changed.**

## References

Jarvis MJ, Belcher M, Vesey C, Hutchison DCS

**Low cost carbon monoxide monitors in smoking assessment**

Thorax 1986; 41: 886-887

Stewart RD, Stewart RS, Stamm W, Seleen RP

### **Rapid estimation of carboxyhemoglobin levels in fire fighters**

JAMA 1976; 235, 390-392

## **Smoke Inhalation/CO Poisoning**

CO poisoning usually occurs as a result of smoke inhalation from fires or exposure to CO from car exhausts or faulty heating systems. CO binds with haemoglobin to form carboxyhemoglobin (COHb) which reduces the capacity of the blood to carry oxygen. Acute poisoning may cause symptoms ranging from headache and breathlessness (at COHb levels of 10% to 30%) through confusion to coma and death (COHb usually greater than 60%). A particularly insidious feature of CO poisoning is the development of neurological problems such as movement disorders (often resembling Parkinson's disease), memory loss and altered personality. Such problems may develop weeks after apparent recovery from acute poisoning.

Chronic exposure to relatively low levels of CO may result in a variety of symptoms including headache, fatigue, poor concentration, dizziness, palpitations, chest pain, visual disturbance, nausea, diarrhoea and abdominal pain.

Chronic CO poisoning as a cause of such symptoms is often missed.

The MicroCO meter allows immediate assessment of patients at risk of CO poisoning who can then be rapidly referred for expert assessment. Prompt treatment with oxygen (in a hyperbaric chamber, if severe) is often life saving. Screening for CO exposure may also reveal the cause of non-specific symptoms relating to low level CO exposure (usually as the result of faulty gas appliances).

## References

Meredith T, Vale A, **Carbon monoxide poisoning**

British Medical Journal, 1988; 296, 77-78

## Cigarette Consumption

The MicroCO meter provides a simple screening test for cigarette consumption for use in anti-smoking clinics and all smoking cessation programs.

Measurement of carboxyhemoglobin has been well validated as an indirect measure of cigarette consumption and is widely used in smoking cessation programs.

Typical values for carboxyhemoglobin and expired CO in smokers and poisoning victims, together with the alarm light status, are given below:

CO(ppm)	%COHb	Cigarette consumption	Indicator
0 – 6	0 – 1	Non smoker	Green
7 –10	1.1 – 1.6	Light smoker	Amber
11 – 72	1.8 – 12	Heavy smoker	Red
>72	>12	Suspected poisoning	Red + alarm

Please note that some urban areas may have high environmental levels of CO. This can cause a rise in exhaled CO of a few ppm above that which is normally present on the breath. In these cases it is possible for a non-smoker to appear at the bottom of the 'light smoker' range (7 – 10 ppm).

## References

### **The relationship between alveolar and blood carbon monoxide concentrations during breath holding**

Jones RH, Ellicott MF, Cadigan JB, Gaensler EA

Journal of Laboratory and Clinical Medicine 1958; 51, 553 – 564

## **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels**

Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369

## **Definition of a reliable threshold value for detecting current smokers by CO measurement**

Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano - Respiratory and TB Unit-General  
Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-ITALY).

## **Correlation between exhaled CO measurements and carboxyhemoglobin percentage in smokers**


Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla; Zanette Antonia;  
Nardini Stefano - Respiratory and TB Unit- General Hospital- Via  
Forlanini, 71- I-31029-Vittorio Veneto (TV- ITALY).

## **Package Contents**


The MicroCO meter is supplied with a carrying case containing the following items:

1. MicroCO meter (Cat No. 36-MC02-STK)
2. 22 mm mouthpiece adapter (Box of 10 Cat No.36-PSA2000)
3. 22 mm reducing connector (Cat No.36-MEC1007)
4. PP3 Battery (Cat No.36-BAT1002)
5. 4 Disposable mouthpieces (Box 250 Cat No.36-PSA1200)
6. Calibration tool (Cat No. 36-MEC1184)


1

2

3

4

5

6


## Operation


For accurate results the CO meter should be used at room temperature. If the instrument has been stored in cool or hot conditions then allow time to reach room temperature prior to use

Install the PP3 battery by sliding open the battery cover, clipping the battery in place and replacing the cover. Insert the mouthpiece adapter into the MicroCO meter and then insert a disposable mouthpiece into the adapter. Turn the unit on by selecting the upper or middle position on the central slide switch and the display will show the software version number:

This will appear momentarily while correction for ambient levels is executed. During this time the unit must not be exposed to raised concentrations of CO. The buzzer will then sound and the display will change to:

Upon hearing the buzzer instructs the subject to inspire fully and hold their breath for 20 seconds. The display will count down from 20 to 0 as an aid to timing the breath holding period. The green indicator will then illuminate and the display will change to:

At this time the subject should seal their lips around the mouthpiece and exhale slowly and fully. The 20 second breath holding period is recommended to allow time for equilibration of alveolar gas.

If, however, the subject is unable to hold their breath for this period, the unit may still be used before 20 seconds have elapsed.

The unit must not be used, however, for one second after the unit is first switched on i.e. before the countdown commences. Expired alveolar gas is then entrapped between sensor and mouthpiece valve. The display reading will rise to a plateau over the course of a several seconds. The final value will be held until the unit is turned off and represents parts per million CO or %COHb depending on the slide switch position.

The lights will come on according to the table on page 4 and the red light will flash and an alarm will sound if the measurement rises above 72ppm (12%COHb).


If this occurs then the possibility of CO poisoning should be investigated as this level of CO is unlikely to be produced by cigarette smoking.

### **Important note:**

**Before repeating a measurement the unit must be turned off, and the mouthpiece and adapter removed for at least 1 minute. This is to allow re-equilibration with ambient air and to dry the surface of the sensor. Visually inspect that all moisture has evaporated from the surface of the sensor before reuse.**

If the unit is switched on again too quickly after use there may be a response to residual expired carbon monoxide from the previous test.

In this case the display will show:

If this is displayed then turn the unit off, remove the mouthpiece adapter, and expose to ambient air for 2 minutes before repeating the test.

**Note:** If this warning appears again after following the above procedure then turn off and leave the sensor exposed to ambient air for a further 3 minutes. If the same message appears again then this indicates possible contamination of the sensor with a solvent.

In this case remove all sources of solvent from the vicinity of the sensor and expose to ambient air for 24 hours before switching on again.

## PC Connection

The MicroCO may be connected to the serial port of a PC running COBRA (CO BReath Analysis) software.

This software allows the measurements to be read by a PC and automatically entered onto a pre-defined report for subsequent printing and filing. It also allows the CO level indicators and the breath holding countdown timer to be configured to individual requirements.


Note: the MicroCO should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 – 'Safety of

Information Technology Equipment including Electrical Business Equipment’.

## Calibration


Calibration will remain stable to within 2% over one month and typically to within 10% over 6 months. CareFusion supplies calibration gas (20ppm CO in air) and recommends that the unit is recalibrated on a 6 monthly basis. See page 35 for calibration accessories.

To carry out the calibration locate the calibration switch on the right hand side of the instrument as shown below.

Push the slide switch to the CO - PPM position and wait for the unit to display zero.

Screw the control valve firmly onto the cylinder and connect the gas supply as shown below:

The plastic tubing supplied with the gas should be pushed firmly over the reducing connector.

Slowly turn the control knob counter-clockwise until the ball in the flow indicator is on the first mark. This will supply the proper amount of

flow to the cell. Apply this flow for 25 seconds and if the meter does not read 20ppm use the calibration tool to push the calibration switch.

**CAUTION: opening the control knob more than indicated will give erroneous reading and could damage the cell.**

The unit will then beep 3 times, store the new calibration value, and display the following:

The gas supply should then be turned off.

If the signal from the CO sensor is too low a new calibration value will not be stored and the display will show:

The most likely cause for this is an expired cell but may also be caused by depressing the calibration switch with no gas applied. Ensure that the concentration of calibration gas is correct (20ppm), that the connections to the gas cylinder are secure, and that the gas cylinder is not empty, and then repeat the calibration procedure.

If the above message is repeated, the CO meter must be returned to CareFusion, for sensor replacement. Sensor life ranges from 2 to 5

years and depends upon both the amount of exposure to CO and other gases, particularly solvents such as alcohol and cleaning fluids.

If the signal from the CO sensor is too high, a new calibration value will not be stored and the following will be displayed when calibration is attempted:

The most likely cause for this is using an incorrect concentration of calibration gas. Ensure that the concentration of calibration gas is correct (20ppm) that the connections to the gas cylinder are secure, and then repeat the calibration procedure.

In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.


### **Important Notes:**

- Only certified calibration gas from a reputable source should be used.
- Ensure that no CO is present on the sensor for 3 minutes before starting the calibration procedure.
- Ensure that the instrument and gas cylinder have stabilized at room temperature before calibrating.
- Do not immerse the unit in the calibration gas.

## **Battery Life**

Battery life is approximately 30 hours of continuous use. When the battery has approximately 1 hour of useful life left the alarm will sound

momentarily after the unit is first switched on and the following message will be displayed:


When the battery has completely expired the above will be displayed continuously and the battery must be replaced.

Note: Please remove the battery if the meter is likely to be unused for more than a month.

### **Internal Battery Expiry**

The MicroCO meter has an internal battery with a life of approximately 10 years. This battery supplies the sensor signal conditioning circuit continuously to ensure instant start up.

When the battery has expired, the following warning message will be displayed:

and the alarm will sound when the unit is first turned on.

When this occurs, the CO meter must be returned to Micro Direct Inc, or an authorized agent, for battery replacement.

## Cleaning

The mouthpiece adapter may be cleaned using a mild detergent solution or replaced if there is risk of contamination. Exposed surfaces of the meter, other than the sensor surface, may be wiped with a clean, damp cloth. It is recommended that this procedure is performed after each use and all used cardboard mouthpieces discarded.

**Important note:** The sensor surface must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol or permanent damage may result.

## Servicing

If your unit requires service or repair please see page 75 for contact details. A full service manual including circuit diagram and parts list is available upon request.

## Consumables / Supporting Products

Cat. No.	Description
36-PSA1200	22mm disposable cardboard mouthpieces (250 per box)
36-PSA2000	22mm mouthpiece adaptor (pack of 10)
36-BAT1002	Alkaline PP3 Battery
36-MEC1184	Calibration Tool
36-MCG020	MicroCan calibration gas. 20 litres of gas in a 1 litre canister. 20ppm Carbon Monoxide, balance air.
36-MGA222	Stainless steel valve for calibration canisters with flow Indicator
36-MEC1007	22mm reducing connector for calibration
36-CAB1000	Interface cable for COBRA software

To place an order for consumables / supporting products, or for general and service enquiries please contact CareFusion or your local dealer.

**For US Customer Care:** Toll Free 1-800-231-2466, phone 1-714-283-2228

## Caution

- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.”

## Symbols

Type B device

## Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992.

## Specifications

Sensor type	Electro-chemical fuel cell
Range	0 - 100 ppm
Resolution	1 ppm
Green indicator light	0 to 6ppm (0 to 1 %COHb)
Amber indicator light	7 to 10ppm (1.1 to 1.6 %COHb)
Red indicator light	11 to 72ppm (1.8 to 12 %COHb)
Flashing red light + alarm	>72ppm (>12 %COHb)
Accuracy	+/-5% of reading or 1ppm whichever is the greater
Sensitivity drift	0.5%/°C
Sensor life	2 to 5 years
Response time	< 15 sec (to 90% of reading)
Hydrogen cross sensitivity	<15%
Operating temperature	15 – 25 °C
Operating pressure	Atmospheric +/- 10%
Pressure coefficient	0.02% signal per mBar
Relative humidity (Non condensing)	15 - 90% continuous (0 - 99% intermittent)
Baseline drift	0ppm (auto-zero)
Long term drift	< 2% signal loss per month
Power source	Single Alkaline 9 volt PP3
Main battery life	30 hours of continuous use equivalent to approximately 2000 tests
Internal battery life	10 years
Weight	180 g (Including battery)
Dimensions	170 x 60 x 26 mm
Display	3 ½ digit LCD
Storage temperature	-20° to +70 °C
Storage humidity	30% to 90%

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## **Référence**

Meredith T, Vale A. **Carbon monoxide poisoning.**

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Jones RH, Ellicott MF, Cadigan JB, Gaensler EA  
Journal of Laboratory and Clinical Medicine 1958; 51, 553 – 564

### **Carbon monoxide in breath in relation to smoking and carboxyhaemoglobin levels**

Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369

### **Definition of a reliable threshold value for detecting current smokers by CO measurement**




Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano - Respiratory and TB Unit-General  
Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-ITALY).

### **Correlation between exhaled CO measurements and carboxyhaemoglobin percentage in smokers**

Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla; Zanette Antonia;  
Nardini Stefano - Respiratory and TB Unit- General Hospital- Via  
Forlanini, 71- I-31029-Vittorio Veneto (TV- ITALY).

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### **Mit Bezug auf:**

#### **“The relationship between alveolar und carbon monoxide concentrations during breath holding”**

Jones RH, Ellicott MF, Cardigan JB, Gaensler EA

Journal of Laboratory and Clinical Medicine 1958; 51, 553-564

#### **“Carbon monoxide in breath in relation to smoking and carboxyhaemoglobin levels”**

Wald NJ, Idle M, Boreham J, Baily A

Thorax 1981; 36, 336-369

#### **“Definition of a reliable threshold value for detecting current smokers by CO measurement”**

Marino Luigi, Latini Roberto, Barbano Gina, Bazzlerla Giorgio, DeLuca

Anita, Nardini Stefano – Respiratory and TB Unit General

Hospital – Via Forlanini, 71-1-31029 – Vittorio Veneto (TV-ITALY).

#### **“Correlation between exhaled CO measurements and carboxyhaemoglobin percentage in smokers”**

Marino Luigi, Latini Roberto, Barbano Gina, Bazzlerla, Zanette Antonia,

Nardini Stefano – Respiratory and TB Unit – General Hospita – Via

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## Customer contact information

### **UK Customers only**

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd  
UK Customer Service & Support  
The Crescent  
Jays Close  
Basingstoke  
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: [micro.uksales@carefusion.com](mailto:micro.uksales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: [micro.ukservice@carefusion.com](mailto:micro.ukservice@carefusion.com)

Technical Support Enquiries:

Telephone: 01256 388551

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

### **International customers only**

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH  
Customer Service & Support International  
Leibnizstrasse 7  
97204 Hoechberg  
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: [micro.internationalsales@carefusion.com](mailto:micro.internationalsales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: [support.admin.eu@carefusion.com](mailto:support.admin.eu@carefusion.com)

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

CareFusion UK 232 Ltd.,  
Quayside  
Chatham Maritime  
Kent ME4 4QY  
U.K.

JDE no. 36-MAN1279  
Drg no. 073-04WW  
Issue 1.0  
January 2010  
© CareFusion 2010